Zuhause TEST

Gesunder Darm

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2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med.(1990), 322: 909–916.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292–296.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am J. Med. 1996; 100:355–415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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INSTRUCTIONS FOR USE

INTENDED USE

The **ZuhauseTEST Healthy Stomach** is a rapid chromatographic immunoassay for the qualitative detection of *H. pylori* antigens in human feces specimens, providing results in 10 minutes. The test utilizes antibodies specific for *H. pylori* antigens to selectively detect *H. pylori* antigens in human feces specimens.

INTRODUCTION

H. pylori is a small, spiral-shaped bacterium that lives in the surface of the stomach and duodenum. It is implicated in the etiology of a variety of gastrointestinal diseases, including duodenal and gastric ulcer, non-ulcer dyspepsia as well as active and chronic gastritis.^{1,2} Both invasive and non-invasive methods are used to diagnose *H. pylori* infection in patients with symptoms of gastrointestinal disease. Specimen-dependent and costly invasive diagnostic methods include gastric or duodenal biopsy followed by urease testing (presumptive), culture, and/or histologic staining.³ A very common approach to the diagnosis of *H. pylori* infection is the serological identification of specific antibodies in infected patients. The main limitation of serology test is the inability to distinguish current and past infections. Antibody may be present in the patient's serum long after eradication of the organisms.⁴ HpSA (*H. pylori* Stool Antigen) testing is gaining popularity for diagnosis of *H. pylori* infection and also for monitoring the efficacy of the treatment of *H. pylori* infection. Studies have found that more than 90% of patients with duodenal ulcer and 80% of patients with gastric ulcer are infected with *H. pylori*.⁵

The **ZuhauseTEST Healthy Stomach** is a qualitative, lateral flow immunoassay for the detection of *H. pylori* antigens in human feces specimens. In this test, the membrane is pre-coated with anti-*H. pylori* antibodies on the test line region of the test. During testing, the specimen reacts with the particle coated with anti-*H. pylori* antibodies. The mixture migrates upward on the membrane by capillary action to react with anti-*H. pylori* antibodies on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

TEST CONTENTS

- 1 test cassette in a sealed pouch
- 1 sample collection applicator with sample dilution buffer
- 1 stool collection paper
- 1 instruction pamphlet

Additionally required:

- 1 timer

STORAGE AND STABILITY

The test cassette can be stored at room temperature or refrigerated (2–30°C). The test cassette is stable through the expiration date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.


TEST PREPARATION

Bring the test cassette and sample dilution buffer at a room temperature (15°C to 30°C) before performing the test.


Please read all the information in this package insert before performing the test.

TEST PERFORMANCE

1. Wash your hands with soap and rinse with clear water. To collect fecal specimens: The stool specimen should be collected in the stool collection paper. Please use the stool collection paper, avoiding contamination of the specimen by taking precautions that the specimen or side of paper containing specimen does not come in contact with any contaminating objects including toilet cleaners.
2. Take the sample collection applicator and open it by screwing the blue cap carefully. Be careful not to spill solution from the tube.


3. Randomly stab the specimen collection applicator into the fecal specimen in at least **3 different sites**. Do not scoop the fecal specimen.
4. Screw on and tighten the cap onto the specimen collection tube, then shake the specimen collection tube vigorously to mix the specimen and the extraction buffer.


5. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the test is performed immediately after opening the foil pouch.

6. Remove the white protective cap at the end of the tube. Using a piece of tissuepaper, break the tip of the tube. Be careful not to spill or spatter the sample mixture from tube.

7. Invert the specimen collection tube and transfer **2 full drops of the extracted specimen** to the specimen well (S) of the test cassette. Avoid trapping air bubbles in the specimen well (S). Then start the timer.

8. Read results at **10 minutes**. Do not read results after 20 minutes.


EVALUATION OF TEST RESULTS

To read the test results simply determine whether a line is present or absent at the control (C) position. It does not matter how strong or weak a Control line is. Staining of the control line (C) serves as an internal progress check. It confirms sufficient sample application and correct test performance.

POSITIVE TEST RESULT


The test result is **positive** if a **colored line appears in the control field (C) and a colored line is visible in the test field (T)**. This result means that *H. pylori* antigen is present in the stool and you should consult a physician.

NOTE: The intensity of the color in the test line region (T) will vary depending on the concentration of *H. pylori* antigen present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive. The final diagnosis must be confirmed by a physician.

NEGATIVE TEST RESULT


The test result is **negative** if a **colored line appears in the control field (C) and no colored line is visible in the test field (T)**. This result means that the presence of the *H. pylori* antigen in feces was not detectable. If there are any other gastrointestinal symptoms, please, visit a physician for further medical investigation.

INVALID TEST RESULT

If there is **no Control line (C) or only a Test line (T)** in the result window, the test did not run correctly and the results are not valid.

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

PERFORMANCE EVALUATION

Sensitivity and Specificity
The **ZuhauseTEST Healthy Stomach** has been evaluated with specimens obtained from a population of symptomatic and asymptomatic individuals. The result shows that the sensitivity of the **ZuhauseTEST Healthy Stomach** is 97.6% and the specificity is 97.9% relative to other rapid test.

Method	Other Rapid Test		Total Results
	Positive	Negative	
Zuhause-TEST Healthy Stomach	Positive	83	2
	Negative	2	93
Total Results		85	95

Relative Sensitivity: 97.6% (95%CI:*91.8%-99.7%)
Relative Specificity: 97.9% (95%CI:*92.6%-99.7%)
Overall accuracy: 97.8% (95%CI:*94.4%-99.4%)
*Confidence Intervals

Intra-Assay
Within-run precision has been determined by using 15 replicates of four specimens: negative, low titer positive, middle titer positive and high titer positive specimens. The specimens were correctly identified >99% of the time.

Inter-Assay
Between-run precision has been determined by 15 independent assays on the same four specimens: negative, low titer positive, middle titer positive and high titer positive specimens. Three different lots of the **ZuhauseTEST Healthy Stomach** have been tested using these specimens. The specimens were correctly identified >99% of the time.

Cross-reactivity
Cross reactivity with following organisms has been studied at 1,0*10⁸ organisms/mL. The following organisms were found negative when tested with the **ZuhauseTEST Healthy Stomach**: *Acinetobacter calcoaceticus*, *Acinetobacter spp*, *Branhamella catarrhalis*, *Candida albicans*, *Chlamydia trachomatis*, *Enterococcus faecium*, *E.coli*, *Enterococcus faecalis*, *Gardnerella vaginalis*, Group A *Streptococcus*, Group B *Streptococcus*, Group C *Streptococcus*, *Hemophilus influenzae*, *Klebsiella pneumonia*, *Neisseria gonorrhoea*, *Neisseria meningitidis*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, Rotavirus, *Salmonella choleraesuis*, *Staphylococcus aureus*, Adenovirus

Interfering Substances

The following potentially Interfering Substances were added to HPG negative and positive specimens. Ascorbic acid: 20 mg/dL; Oxalic acid: 60 mg/dL; Bilirubin: 100 mg/dL; Uric acid: 60 mg/dL; Aspirin: 20 mg/dL; Urea: 2000 mg/dL; Glucose: 2000 mg/dL; Caffeine: 40 mg/dL; Albumin: 2000 mg/dL

WARNINGS AND IMPORTANT INFORMATION


- For self-testing *in vitro* diagnostic use only. The test is for *in vitro* diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use the after expiration date.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.
- Store in a dry place at 2–30°C (36–86°F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- If you use a container to collect the stool samples, it must be clean.
- Follow the indicated time strictly.
- Use the test only once. Do not dismantle and touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- The test should be used for the detection of *H. pylori* antigens in feces specimens only. Neither the quantitative value nor the rate of increase in *H. pylori* antigens concentration can be determined by this qualitative test.
- The **ZuhauseTEST Healthy Stomach** will only indicate the presence of *H. pylori* in the specimen and should not be used as the sole criteria for *H. pylori* to be etiologic agent for peptic or duodenal ulcer.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of *H. pylori* infection.
- Following certain antibiotic treatments, the concentration of *H. pylori* antigens may decrease to the concentration below the minimum detection level of the test. Therefore, diagnosis should be made with caution during antibiotic treatment.

LITERATURE


1. Marshall, BJ, McGeechie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439–444.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med.(1990), 322: 909–916.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292–296.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am J. Med. 1996; 100:355–415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

EXTRA INFORMATION

1. **How does the ZuhauseTEST Healthy Stomach work?** *H. pylori* is a small, spiral-shaped bacterium that lives in the surface of the stomach and duodenum. The **ZuhauseTEST Healthy Stomach** detects specifically the antigens in feces to ascertain the presence of the bacterium.
2. **When should the test be used?** The test can be performed anytime of the day. The test can be performed in case of repeated stomach and intestinal troubles (GERD, gastritis etc.).
3. **Can the result be incorrect?** The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if **ZuhauseTEST Healthy Stomach** gets wet before performing the test or if the quantity of feces dispensed in the sample well is too much or not sufficient, or if the number of extracted specimens drops are less than 2 or more than 3. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A

CONTENTS OF THE KIT
* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the “content” section on the external box.

- Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
- Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.
- The desiccant packet must not be used. Dispose of it with household waste without opening it.

- TEST PROCEDURE**
- The test must be used immediately after opening the sealed aluminium pouch. If the device has been stored at refrigerated temperature (4 °C), allow it to reach room temperature for at least 30 minutes before use.
- Wash hands with soap and warm water, rinse with clean water and allow to dry. Note: The use of warm water facilitates capillary blood collection as it induces vasodilation.
 - Lay out the material needed as follows: open the aluminium bag, take out only the test cassette and dispose of the desiccant packet. Open the plastic packet containing the pipette. Open the vial by unscrewing the white cap, making sure it does not fall to the ground. –**FIG. A**
 - Carefully rotate the protective cap of the sterile lancet 360° without pulling it, then extract and discard the cap. –**FIG. B**
 - Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended). It is important that the massage is done from the palm of the hand to the fingertip, to improve blood flow. Press the open end of the lancet (the side the cap has been extracted from), against the fingertip –**FIG. C**. The tip of the lancet automatically retracts after use. If the lancet does not work properly, discard it and use the second one supplied. If the second one is not required, it can be disposed of without special precautions.
 - Holding the hand downwards, massage the finger until a large drop of blood forms. It is important to massage from the palm of the hand to the fingertip to improve blood flow. –**FIG. D**
 - Take the pipette without pressing the bulb. Two sampling methods are suggested:
–**FIG. E1**: hold the pipette horizontally without pressing the bulb place it in contact with the drop of blood, it will enter the pipette by capillarity. Move the pipette away when the black line is reached. If there is not enough blood, continue to massage the finger until the black line is reached.
–**FIG. E2**: place the pipette on a clean, flat surface with the tip protruding from the surface, then place the drop of blood in contact with the pipette, it will enter by capillarity. If blood is not sufficient, continue massaging the finger until the blood has reached the black line.
Avoid, as far as possible, to continuously move the tip of the pipette away from the finger in order to prevent the formation of air bubbles.
 - Deposit the blood collected with the pipette into the previously opened vial. Press the pipette bulb 2 or 3 times to ensure that all the blood has moved into the diluent solution. Replace the screw cap onto the vial and mix well for at least 10 seconds. –**FIG.F**
 - Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). –**FIG.G**. Deposit 3 drops into the well shown on the cassette (S). –**FIG.H**
 - Wait 5 minutes and then read the results as indicated.

RESULTS INTERPRETATION
READ THE RESULT AT EXACTLY 5 MINUTES.
DO NOT READ THE RESULT AFTER 6 MINUTES
The result must be interpreted independently of the colour intensity of the bands that may appear in the reading window.

NEGATIVE, C-REACTIVE PROTEIN BELOW 10 MG/L
A coloured band appears only under the C (Control) sign, as shown in the figure on the right. This result indicates that the concentration of CRP is below 10 mg/L and suggests the absence of infection and inflammation.

POSITIVE, C-REACTIVE PROTEIN BETWEEN 10 AND 40 MG/L
Two coloured bands appear in the reading window, one by the C (Control) sign, the other by the T1 sign, as shown in the figure on the right. It is recommended to consult a physician, as CRP values between 10 and 40 mg/L may be associated with a viral infection or the onset of a bacterial infection. An inflammation caused (for example) by physical trauma can also lead to this kind of result.

POSITIVE, C-REACTIVE PROTEIN BETWEEN 40 AND 80 MG/L
In the reading window, three coloured bands appear, one by the C (Control) sign and two others by the T1 and T2 signs, as shown in the figure on the right. Consult a physician, as this result may be associated with a viral or bacterial infection or physical trauma.

POSITIVE, C-REACTIVE PROTEIN OVER 80 MG/L
Four coloured bands appear in the reading window, one by the C (Control) sign and three others by the T1, T2 and T3 signs, as shown in the figure on the right. Consult a doctor, as this result may be associated with a bacterial infection or severe inflammation.
NOTE: The colour intensity in the control and test line regions may vary (as represented in the images). The intensity of the bands is not relevant; only their presence must be considered.

NON VALID RESULT
The control band does not appear (no lines appear or one or more T bands appear without the C band).
Insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the control line. Review the procedure and repeat the test with a new device and a new blood sample.

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EN Iron Deficiency

INTRODUCTION

The **ZuhauseTEST Iron Deficiency** can be used as a tool for the detection of iron deficiency. Anemia due to iron depletion is widely held in children and women of all ages but mainly by women who still have their periods (at least 20% suffer from iron deficiency). Main signs are paleness, feeling tired, headaches, faster heartbeat, or shortness of breath during exercise. They may appear gradually and could go unnoticed. Iron deficiency occurs when blood does not contain enough red blood cells and thus low levels of hemoglobin, which is the major protein involved in oxygen transport in whole body. An important component of hemoglobin is iron.

Depletion of iron, which can happen during pregnancy, growth, in case of insufficient iron intake, inadequate absorption or blood loss (period, abnormal bleedings, ulcers, etc.) has tremendous effects on health. Low ferritin may also indicate hypothyroidism, vitamin C deficiency or celiac disease. Low ferritin levels are seen in some patients with restless legs syndrome, not necessarily related to anemia, but perhaps due to low iron stores just prior to anemia.^{1,2}

The **ZuhauseTEST Iron Deficiency** is a rapid chromatographic immunoassay for the qualitative detection of ferritin in human fingertip blood at a cut-off concentration of 30 ng/mL. The membrane is pre-coated with anti-ferritin polyclonal antibody on the test line region. The gold layer is pre-coated with anti-ferritin monoclonal antibody and Rabbit IgG. During testing, the specimen reacts with the particle coated with anti-ferritin monoclonal antibody. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-ferritin polyclonal antibody on the membrane and generate a colored line. The line in test line region (T) appears, if the ferritin level exceeds the cut-off level of 30 ng/mL. If the ferritin concentration is less than 30 ng/mL, the test line does not appear. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

Analytical Sensitivity: With the **ZuhauseTEST Iron Deficiency** ferritin levels in human fingertip blood as low as 30 ng/mL can be detected.

Cross-Reactivity: An evaluation was performed to determine the cross reactivity and interferences of **ZuhauseTEST Iron Deficiency**. There is no cross reactivity with HAMA, RF, Human serum albumin, human AFP, Ferric Chloride, human transferrin and human hemoglobin.

STORAGE AND STABILITY

Store the test in the sealed foil pouch at room temperature or refrigerated (2–30 °C). The test is usable until the expiration date printed on the sealed foil pouch. **DO NOT FREEZE**. Do not use after the expiration date.

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- For self-testing, *in vitro* diagnostic use only. The test is for external use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Keep out of the reach of children.
- All test components can be disposed of in household waste.
- Store in a dry place at 2–30 °C, avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- This test kit is intended to be used as a preliminary test only and repeatedly abnormal results should be discussed with doctor or medical professional.

1 Alcohol swab
Jiangsu Sunclan Medical Co., Ltd.
No. 11 Fenghuang south Road, Hutang Town, Wujin District 213162 Changzhou City, Jiangsu Province, P.R. China.

oder
LIGHTS MEDICAL MANUFACTURE CO., LTD.
No. 19 Quanda Road, Wujing Development Area, Tianjin, China

CE 0197

oder
Baoying Country Fukang Medical Appliance Co., Ltd.
Guangyang road, Huangcheng town industrial area, Baoying County, Yangzhou, Jiangsu, 225800, China

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PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Cangxing Street, Cangxian Community, Yuhang District, Hangzhou City, 311212 Zhejiang, China

CE 0197

1 Compresse d'alcool
Jiangsu Sunclan Medical Co., Ltd.
No. 11 Fenghuang south Road, Hutang Town, Wujin District 213162 Changzhou City, Jiangsu Province, P.R. China.

ou
LIGHTS MEDICAL MANUFACTURE CO., LTD.
No. 19 Quanda Road, Wujing Development Area, Tianjin, China

CE 0197

ou
Baoying Country Fukang Medical Appliance Co., Ltd.
Guangyang road, Huangcheng town industrial area, Baoying County, Yangzhou, Jiangsu, 225800, China

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
EN Wash your hands with soap and rinse them with clear, warm water. Bring the foil pouch containing the test to room temperature before opening it. Open the pouch, take out the test cassette and the pipette, and place them on a clean and flat surface.

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
EN Carefully twist off the cap of the lancing device and dispose of it.

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EN Use the provided alcohol pad to clean the fingertip of the middle or ring finger as the puncture site. Allow to air dry.

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EN Press the lancing device against the finger on the side where the cap was removed. After use, the tip will automatically and safely retract.

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EN Dispense the entire blood sample into the sample well (S) of the cassette by squeezing the bulb of the pipette.

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EN Read the results exactly after 5 minutes. The result must not be interpreted after 10 minutes.

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EN Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and add 1 drop of buffer into the Buffer well (S) of the cassette.

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LITERATUR

1. Merck Manual of Diagnosis and Therapy, Thyroid gland disorders.
2. The American Heritage Dictionary of the English Language, Fourth Edition. Houghton Mifflin Company. 2006. ISBN 0-395-82517-2.
3. Sacher R, Richard A. McPherson (2000). Widmann's Clinical Interpretation of Laboratory Tests, 11th ed. F.A. Davis Company. ISBN 0-8036-0270-7.
4. So, M; MacIsaac, RJ; Grossmann M (August 2012). „Hypothyroidism“. Australian Family Physician 41 (8): 556–62.
5. Surkset. al. JAMA 291:228, 2004 Daniel, GH, Martin, JB, Neuroendocrine Regulation and Diseases of the Anterior Pituitary and Hypothalamus in Wilson, JD, Braunwald, E., Isselbacher, KJ, et. al., Harrison's Principles of Internal Medicine, 12th Edition, McGraw-Hill, Inc., New York, NY, 1991, p. 1666)

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EC REP
MedNet EC-REP GmbH
Borkstrasse 10, 48163 Münster
Germany

CE 0123

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PROMISEMED HANGZHOU
MEDITECH CO., LTD.
No. 1388 Cangjiang Street,
Cangjian Community
Yuhang District, Hangzhou City,
311127 Zhejiang, China

CE 0197

oder

Ningbo Medsun Medical Co., Ltd.
No. 298 Huangjiu Road,
Jiangbei 315031 Ningbo,
People's Republic of China

CE 0123

oder

Shandong Lianfa Medical Plastic
Products Co., Ltd.
No. 1 Shuangshan Sanjian Road,
250200, Zhangqiu City, Jinan,
Shandong,
PEOPLE'S REPUBLIC OF CHINA

CE 0123

Alkoholtupfer

Baoying County Fukang Medical
Appliance Co., Ltd.
Guangyang road, Huangcheng
town industrial area, Baoying
County, Yangzhou, Jiangsu,
225800, China

CE 0197