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

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- * Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
- * Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.
- * The desiccant packet must not be used. Dispose of it with household waste without opening it.

TEST PROCEDURE

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
Note: The use of warm water facilitates capillary blood collection as it induces vasodilation.
- 2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag. Open the plastic packet containing the pipette. –**FIG. A**
- 3) Carefully rotate the protective cap of the sterile lancet 360° without pulling it. Extract and discard the released cap. –**FIG. B**
- 4) Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended). It is important that the massage is done from the palm of the hand to the phalanx, to improve blood flow. Press the open end of the lancet (the side the cap has been extracted from), against the fingertip –**FIG. C** The tip of the lancet automatically retracts after use. If the lancet does not work properly, discard it and use the second one supplied. If the second one is not required, it can be disposed of without special precautions.
- 5) Holding the hand down, massage the finger until a large drop of blood forms. It is important to massage from the palm of the hand to the phalanx to improve blood flow. –**FIG. D**
- 6) Take the pipette without pressing the bulb. Two sampling methods are suggested:
–**FIG. E1:** hold the pipette horizontally without pressing the bulb place it in contact with the drop of blood, it will enter the pipette by capillarity. Move the pipette away when the black line is reached. If there is not enough blood, continue to massage the finger until the black line is reached.
–**FIG. E2:** place the pipette on a clean, flat surface with the tip protruding from the shelf, then place the drop of blood in contact with the pipette, it will enter by capillarity. If blood is not sufficient, continue massaging the finger until the blood has reached the black line.
Avoid, as far as possible, to continuously move the tip of the pipette away from the finger in order to prevent the formation of air bubbles.
- 7) Place the blood collected with the pipette into the well indicated on the cassette (S) by pressing the pipette bulb. –**FIG. F**
- 8) Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). Deposit 2 drops into the well indicated on the cassette (S), wait for 5 seconds between the first and second drop. –**FIG. G**
- 9) Wait 5 minutes and read the results referring to the next section for result interpretation.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES

The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line.

This result means that the value of Prostate-Specific Antigen (PSA) in the blood is higher than normal, so you must consult a physician.

NEGATIVE RESULT

A coloured band appears only under the C (control) sign.

This means that the value of Prostate-Specific Antigen (PSA) in the blood is normal.

A negative result does not exclude the possibility of prostatic pathologies.


NON VALID RESULT

No bands appear or there is a line only under the T (test) sign and not under the C (control) sign. In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test with a new blood sample.

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EN Vitamin D

INTRODUCTION

Vitamin D refers to a group of fat-soluble secosteroids responsible for increasing intestinal absorption of calcium, iron, magnesium, phosphate and zinc. In humans, the most important compounds in this group are vitamin D3 and vitamin D2.¹ Vitamin D3 is naturally produced in the human skin through the exposure to ultraviolet light and Vitamin D2 is mainly obtained from foods. Vitamin D is transported to the liver where it is metabolized to 25-hydroxy Vitamin D. In medicine, a 25-hydroxy Vitamin D blood test is used to determine Vitamin D concentration in the body. The blood concentration of 25-hydroxy Vitamin D (including D2 and D3) is considered the best indicator of Vitamin D status. Vitamin D deficiency is now recognized as a global epidemic.² Virtually every cell in our body has Receptors for Vitamin D, meaning that they all require "Sufficient" Level of Vitamin D for adequate functioning. The health risks associated with Vitamin D deficiency are far more severe than previously thought. Vitamin deficiency has been linked to various serious diseases: Osteoporosis, Osteomalacia, Multiple Sclerosis, Cardiovascular Diseases, Pregnancy Complications, Diabetes, Depression, Strokes, Autoimmune Diseases, Flu, Different Cancers, Infectious Diseases, Alzheimer, Obesity and Higher Mortality etc.³

The **ZuhauseTEST Vitamin D** is a rapid chromatographic immunoassay for the semi-quantitative detection of 25-hydroxyvitamin D (25 (OH) D) in human fingerstick whole blood. This assay provides a preliminary diagnostic test result and can be used to screening for Vitamin D deficiency. The **ZuhauseTEST Vitamin D** is based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip. During testing, 25 (OH) D present in the specimen will compete with 25 (OH) D on the test line for limited amount of anti-25 OH Vitamin D antibodies in the conjugate. The higher concentration of 25 (OH) D in the specimen, the lighter would be the T line. The result will be read according to the Color card provided with the kit. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and that the test has been conducted correctly.

The **ZuhauseTEST Vitamin D** is a rapid chromatographic immunoassay for the semi-quantitative detection of 25-hydroxyvitamin D (25 (OH) D) in human fingerstick whole blood. This assay provides a preliminary diagnostic test result and can be used to screening for Vitamin D deficiency. The **ZuhauseTEST Vitamin D** is based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip. During testing, 25 (OH) D present in the specimen will compete with 25 (OH) D on the test line for limited amount of anti-25 OH Vitamin D antibodies in the conjugate. The higher concentration of 25 (OH) D in the specimen, the lighter would be the T line. The result will be read according to the Color card provided with the kit. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and that the test has been conducted correctly.

The **ZuhauseTEST Vitamin D** is a rapid chromatographic immunoassay for the semi-quantitative detection of 25-hydroxyvitamin D (25 (OH) D) in human fingerstick whole blood. This assay provides a preliminary diagnostic test result and can be used to screening for Vitamin D deficiency. The **ZuhauseTEST Vitamin D** is based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip. During testing, 25 (OH) D present in the specimen will compete with 25 (OH) D on the test line for limited amount of anti-25 OH Vitamin D antibodies in the conjugate. The higher concentration of 25 (OH) D in the specimen, the lighter would be the T line. The result will be read according to the Color card provided with the kit. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and that the test has been conducted correctly.

The **ZuhauseTEST Vitamin D** is a rapid chromatographic immunoassay for the semi-quantitative detection of 25-hydroxyvitamin D (25 (OH) D) in human fingerstick whole blood. This assay provides a preliminary diagnostic test result and can be used to screening for Vitamin D deficiency. The **ZuhauseTEST Vitamin D** is based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip. During testing, 25 (OH) D present in the specimen will compete with 25 (OH) D on the test line for limited amount of anti-25 OH Vitamin D antibodies in the conjugate. The higher concentration of 25 (OH) D in the specimen, the lighter would be the T line. The result will be read according to the Color card provided with the kit. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and that the test has been conducted correctly.

- For self-testing *in vitro* diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- The result of this test is an indication of the vitamin D status. In case of repeatedly abnormal results or symptoms, a doctor must be consulted.
- Follow the indicated time strictly, otherwise the results may be inaccurate.
- Use the test only once. Do not dismantle and touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.
- The **ZuhauseTEST Vitamin D** provides only a semi-quantitative analytical result. A secondary analytical method must be used to obtain a confirmed result.
- It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen may cause erroneous results.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- Other clinically available tests are required if questionable results are obtained.

TEST CONTENTS

- 1 test cassette [25(OH)D] with 1 pipette in a foil pouch
- 1 tube with sample dilution buffer
- 2 automatic lancet devices

PROMISED MEDITECH CO., LTD.

No. 1388 Gangping Street, Cangshan Community, Yuhang District, Hangzhou City, 311121 Zhejiang, China

- 1 alcohol swab

Baoying County Fukang Medical Appliance Co., Ltd.

Guangyang road, Huangcheng town industrial area, Baoying County, Yangzhou, Jiangsu, 225800, China

- 1 Instruction for use

- 1 color card

Additionally required:

- 1 timer

TEST PREPARATION

Bring the test cassette and sample dilution buffer to room temperature before starting the test. Wash your hands thoroughly with soap and rinse them with clear, warm water.

TEST PERFORMANCE

Read the instructions for use **completely** before performing the test. A **step-by-step instruction** is provided on the next page and describes the test procedure.

PERFORMANCE EVALUATION

Accuracy
A clinical evaluation was conducted comparing the results obtained using the **ZuhauseTEST Vitamin D** with another Vitamin D rapid test. The in-house clinical trial included 90 whole blood specimens. The results demonstrated with an overall accuracy of 94.4%.

Method		Reference Method (Vitamin D Rapid Test)			Total Results
Zuhause-TEST Vitamin D	Results	Deficient	Insufficient	Sufficient	
	Deficient	4	3	0	7
	Insufficient	0	53	2	55
	Sufficient	0	0	28	28
Total Results		4	56	30	90
Accuracy		> 99.9%	94.6 %	93.3 %	94.4 %

STORAGE AND STABILITY

Store the test in the sealed pouch at 2-30 °C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

WARNINGS AND IMPORTANT INFORMATION

- For self-testing *in vitro* diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- The result of this test is an indication of the vitamin D status. In case of repeatedly abnormal results or symptoms, a doctor must be consulted.
- Follow the indicated time strictly, otherwise the results may be inaccurate.
- Use the test only once. Do not dismantle and touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.
- The **ZuhauseTEST Vitamin D** provides only a semi-quantitative analytical result. A secondary analytical method must be used to obtain a confirmed result.
- It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen may cause erroneous results.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- Other clinically available tests are required if questionable results are obtained.

Rapid test for the detection of vitamin D deficiency (whole blood test by finger puncture)

For self-testing

BIOLOGICAL REFERENCE RANGES AND LITERATURE

25-OH Vitamin D Level	Reference Range (ng/mL)	Reference Range (nmol/L)
Deficient	0-10	0-25
Insufficient	10-30	25-75
Sufficient	30-100	75-250

- Holick MF (March 2006). High prevalence of vitamin D inadequacy and implications for health.Mayo Clinic Proceedings. 81 (3): 353–73.

- Eriksen EF, Glerup H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. Biogerontology. 3 (1-2): 73–7.

- Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. Alternative Medicine Review.10 (2): 94–111.

EXTRA INFORMATION

- How does the Vitamin D test work?** In medicine, a 25-hydroxy Vitamin D is the main storage form of vitamin D in the body. Therefore, the overall status of vitamin D can be determined by detecting the content of 25-hydroxy Vitamin D. 25-hydroxy Vitamin D level less than 30 ng/mL in case of a positive result, indicates Vitamin D Deficiency or Insufficiency. Vitamin D supplements can be recommended in these cases.

- When should the test be used?** The clinical application of 25-hydroxy Vitamin D is mainly for diagnosis, treatment and monitoring of rickets (children), osteomalacia, postmenopausal osteoporosis and renal osteopathy. Vitamin D deficiency is also associated with many other diseases, including cancer, cardiovascular disease, autoimmune diseases, diabetes and depression. Monitor your vitamin D levels to determine whether to take vitamin D supplements. The Vitamin D Rapid Test can be used any time of the day.

- Can the result be incorrect?** The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the **ZuhauseTEST Vitamin D** gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3. The pipette included in the foil pouch allows for checking the correct amount of collected blood. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

- How to interpret the test if the color and the intensity of the lines are different?** Please refer to the illustration and compare the T line intensity with "Vitamin D Color card" provided with the kit.

- If I read the result after 20 minutes, will the result be reliable?** No. The result should be read at 10 minutes after adding the buffer. The result is unreliable after 20 minutes.

- What do I have to do if the result is deficient or insufficient?** If the result is deficient or insufficient, it means that the Vitamin D level in blood is less than 30 ng/mL and that you should consult a physician to show the test result. The physician will decide whether additional analysis should be performed.

- What do I have to do if the result is sufficient?** If the result is sufficient, it means that the Vitamin D level is higher than or equal to 30 ng/mL and is within the normal range. A case of Vitamin D toxicity (hypercalcemia), though rare, but cannot be excluded based on such test results. However, if the symptoms persist, it is recommended to consult a physician.

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PROMISEMED HANGZHOU MEDITECH CO., LTD.

No. 1388 Cangxing Street, Cangan Community, Yuhang District, Hangzhou City, 311121 Zhejiang, China

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EN The test result indicates a Vitamin D deficiency if a light to dark red line appears in the control field (C) and a T-line is visible in the test field (T), which is as dark or darker than the 10 ng/mL line on the test card (upper T-line). For a definitive diagnosis, you should consult your doctor.

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EN The test result indicates an insufficient amount of Vitamin D if a light to dark red line appears in the control field (C) and a T-line is visible in the test field (T), which is darker than the 30 ng/mL line (lower T-line) and lighter than the 10 ng/mL line on the test card (upper T-line). The intensity of the T-line corresponds approximately to the „Insufficient“ line depicted on the test card (middle T-line). For a definitive diagnosis, you should consult your doctor.

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EN The test result indicates a sufficient amount of Vitamin D if a light to dark red line appears in the control field (C) and a T-line is visible in the test field (T), which is as intense or lighter than the 30 ng/mL line on the test card (lower T-line). If no T-line appears, this may indicate an excess of Vitamin D. If symptoms suggestive of a deficiency persist, we recommend consulting a doctor.

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EN If no control line (C) appears or only a test line (T) is visible, the test has not been performed correctly and is invalid. Please check whether you have followed all the instructions carefully. Conduct a new test using a fresh blood sample. If the issue persists, discontinue using the test kit immediately and contact your local distributor.

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EN Open the pouch, remove the test cassette as well as the pipette and place it on a clean and level surface. Run the test within one hour and best results will be obtained if the test is performed immediately after opening the foil pouch.

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EN Carefully twist off the cap of the lancing device and dispose of it.

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EN Use the provided Alcohol pad to clean the fingertip of the middle or ring finger as the puncture site. Allow to air dry.

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EN Press the lancing device against the finger on the side where the cap was removed. After use, the tip will automatically and safely retract.

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- Holick MF (March 2006). High prevalence of vitamin D inadequacy and implications for health. Mayo Clinic Proceedings. 81 (3): 353–73.
- Eriksen EF, Glerup H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. Biogerontology. 3 (1-2): 73–7.
- Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. Alternative Medicine Review.10 (2): 94–111.

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EN Hold your hand slightly downward and massage it, without touching the puncture site, towards the fingertip to obtain a drop of blood.

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EN Bring the tip of the pipette into contact with the blood without touching the finger. Fill the capillary up to the indicated line with blood. You may massage again your finger to obtain more blood if the blood does not reach the indicated line. Avoid of air bubbles.

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EN Use the pipette to transfer the **entire blood** sample into the specimen well (S) of the cassette.

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EN Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and add **2 drops of buffer** into the Buffer well (B) of the cassette.

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EN Read results **exactly after 10 minutes**. Compare the T line intensity with “Vitamin D Color card” provided with the kit to get the Vitamin D level in your blood. Do not interpret the result after 20 minutes.

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