Erkältungs-Check

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5-in-1 Atemwegs-Check

Seiten 2-3

COVID-19 / Influenza A+B / RSV / Adenoviren

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Seiten 4-6

ZuhauseTEST Streptokokken A

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COVID-19/FLU A+B/ RSV/Adenovirus Antigen Rapid Test

INTENDED USE

The test is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2, Influenza A virus, Influenza B virus, Respiratory Syncytial Virus (RSV) and Adenovirus antigens in human nasal swab specimens as an aid in the diagnosis of SARS-CoV-2, Influenza A, Influenza B, RSV and Adenovirus infections. The test is intended for individuals who are suspected of being infected with SARS-CoV-2, Influenza A/B, RSV and/or Adenovirus. The test is not automated. For self-testing use. For *in vitro* diagnostic use only.

SUMMARY

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected with the novel coronavirus are the main source of infection; asymptomatic infected people can also be a source of infection. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.¹ Influenza (commonly known as 'flu') is a highly contagious, acute viral infection of the respiratory tract. Virus transmission occurs when a susceptible individual comes into contact with aerosols or respiratory fomites from an infected individual.² Influenza outbreaks occur each year during the fall and winter months. Type A viruses are typically more prevalent than type B viruses and are associated with most serious influenza epidemics, while type B infections are usually milder. Respiratory Syncytial Virus (RSV), which causes infection of the lungs and breathing passages, is a major cause of respiratory illness in young children. In adults, it may only produce symptoms of a common cold, such as a stuffy or runny nose, sore throat, mild headache, cough, fever, and a general feeling of being ill. Most children with RSV infection, both those who were hospitalized and those who were treated as outpatients, had no existing medical conditions or characteristics that significantly identified them as being at greater risk for severe RSV disease, except for being under 2 years of age.³ Human Adenoviruses comprise an important group of etiologic agents that are responsible for various diseases in adults and children, such as respiratory, ocular, gastroenteric, and urinary infections. In immunocompromised and organ-transplanted individuals, these agents can cause generalized infections.⁴

HOW DOES IT WORK


The test is a qualitative, lateral flow immunoassay for the detection of SARS-CoV-2, Influenza A and Influenza B viruses, RSV, and Adenovirus antigens in human nasal swab specimen. When the extracted specimen is added to the specimen well of the test cassette, the SARS-CoV-2, Influenza A virus, Influenza B virus, RSV, or Adenovirus antigens, if present in the specimen, react specifically with the virus antibodies coated onto the particles, forming a mixture. The mixture migrates up the membrane and reacts with the virus antibodies immobilized on the membrane, resulting in the generation of one or two colored lines in the test line region (T/A/B) of each test window. The presence of the colored line(s) in the test line region(s) indicates a positive result, while its absence indicates a negative result. To serve as a built-in quality control, a colored line will always appear in the control line region(C), indicating that the test cassette functions well.

REAGENTS

The test contains Mouse IgG coated colloidal gold particle, Rabbit IgG coated colloidal gold particle, SARS-CoV-2 antibody coated latex particle, Influenza A antibody coated colloidal gold particle, Influenza B antibody coated colloidal gold particle, RSV antibody coated colloidal gold particle, Adenovirus antibody coated colloidal gold particle, and contains Goat anti-Mouse IgG antibody, Goat anti-Rabbit IgG antibody, SARS-CoV-2 antibody, Influenza A antibody, Influenza B antibody, RSV antibody, Adenovirus antibody coated on the membrane.

KIT COMPONENTS

Materials Provided			
Components	Kit size	1 T/klt	5 T/klt
	Test cassette(s)	1	5
	Package insert	1	1
	Sterile swab(s)	1	5
	Buffer 0.55mL (0.02% Proclin 300)	1	5
Tube holder		On the box	

Materials Required But Not Provided

- Timer

SPECIMEN COLLECTION


Please wash and dry your hands thoroughly. Remove the sterile swab from the pouch. Do not touch the soft tip of the swab. Nasal swab specimen collection:

1. Insert the swab into one nostril (Approx. 2 cm deep).
2. Slowly rotate the swab for 5-10 times, rubbing it against the nasal wall.

Note: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain. When the nasal mucosa is damaged or bleeding, nasal swab collection is not recommended. If you are swabbing others, please wear a face mask. With children, you may not need to insert the swab so deep into the nostril. For very young children, you may need another person to steady the child's head while swabbing. 3. Remove the swab and repeat the collection in the other nostril using the same swab, avoiding excessive or highly viscous nasal discharge. **Note:** Test should be performed immediately after the specimen has been collected.

HOW TO DO THE TEST

Allow the test cassette, the specimen and buffer to equilibrate to room temperature (15-30°C) prior to testing.


Remove the cover of the Buffer tube and insert the tube into the dotted circle on the packaging box or into the tube holder.


"Test it quickly"

self-test for use at home

READING THE RESULTS

STORAGE AND STABILITY

1. Store the kit at 2-30 °C. DO NOT FREEZE.
2. Keep the kit away from sunlight, moisture and heat.
3. Do not use the kit beyond the expiration date.
4. Please use the test cassette within **one hour** after removing it from the foil pouch.

WARNINGS AND PRECAUTIONS

Please read all the information in this package insert before performing the test.

1. This kit is for *in vitro* diagnostic use only.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Do not drink the buffer in the kit. Carefully handle the buffer and avoid it contacting skin or eyes. Rinse with plenty of running water immediately if contacts.
4. Do not use the kit if any of the components or packaging is damaged.
5. **Use the kit only once** and follow the test procedures strictly. Do not dismantle the test cassette or touch the test window of the test cassette.
6. Keep the kit out of the reach of children. Testing for children should be conducted by an adult.
7. Handle all specimens as if they contain infectious agents. The used tests and specimens may be infectious and should be discarded according to local regulations.
8. Wash hands thoroughly before and after handling.
9. Please ensure that an appropriate amount of specimen is used for testing. Too much or too little specimen may lead to deviation of results.
10. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.
11. The components provided in the kit are approved for use with this test. Do not use components from any other commercial kit.
12. This test provides only a preliminary result and should not be considered a substitute for medical advice. Always seek guidance from a qualified healthcare provider before making any health-related decisions.

LIMITATIONS

1. The **SPECIMEN COLLECTION, HOW TO DO THE TEST and READING THE RESULTS** must be followed closely. Failure to follow the procedures may lead to inaccurate test results.
2. Children aged 3-14 years should be tested by an adult. Do not conduct the test if you do not feel confident testing a child. If the child experiences pain, stop the test immediately. This kit is not suitable for children under 3 years old.
3. If the test result is negative and clinical symptoms persist, it is probably because the amount of the virus in very early infection stage is too small to be detected, and it is recommended to do the test again with a new kit or with a molecular diagnostic device to rule out infection.
4. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic technique should be considered to rule out infection in these individuals.
5. Positive results of COVID-19 may be due to infection with non-SARS-CoV-2 coronavirus strains or other interfering factors.
6. A positive result does not preclude an underlying co-infection with another pathogen, therefore the possibility of an underlying bacterial infection should be considered.
7. Excess blood or mucus on the swab specimen may interfere with test performance, and even may yield a false positive result.
8. A false negative result may be obtained if the concentration of the viruses present in the specimen is not adequate or below the detectable level of the kit, or if you fail to follow these procedures such as improper specimen collection or testing.
9. The kit will only indicate the presence or the absence of SARS-CoV-2, Influenza A, Influenza B, RSV and Adenovirus antigens in the nasal swab specimen. Neither the quantitative value nor the rate of increase in the concentration of these viruses can be determined by this qualitative test.

QUESTIONS & ANSWERS

1. How does the COVID-19/FLU A+B/RSV/Adenovirus Antigen Rapid Test work?
The kit can qualitatively detect SARS-CoV-2, Influenza A, Influenza B, RSV and Adenovirus antigens in nasal swab specimens through the specific antibodies it contains. A positive result indicates that SARS-CoV-2, Influenza A, Influenza B, RSV and/or Adenovirus is/are probably present in the specimen.

2. When should the kit be used?

You can do this test as an early diagnosis of SARS-CoV-2/Influenza A+B/RSV/Adenovirus infection when you (or the patient) have (has) the symptoms such as fever or chills, cough, fatigue, sore throat, headache, difficulty breathing, runny or stuffy nose, muscle and joint pain, wheezing, sneezing, decrease in appetite, new loss of taste or smell, nausea or vomiting, diarrhea, severe malaise (feeling unwell), bronchospasm. You can also use this test kit when you suspect that you (or the patient) may have been infected even you (or the patient) don't (doesn't) have any symptoms since there exists the chance of asymptomatic case.

3. Can the result be incorrect?

The result may be incorrect if inadequate or improper specimen is collected or the COVID-19/FLU A+B/RSV/Adenovirus Antigen Rapid Test gets wet before being used, or if the number of drops of solution added is less than 3 or more than 4. Besides, due to the immunological principle involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

4. How do I interpret the test results if the color and the intensity of the lines are different?

The color and the intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible. It should be considered as positive whatever the color or intensity of the test line is.

5. My test result is negative. Does that mean I'm not infected?

A negative result means that you (or the patient) are (is) probably not infected or that the viral load is too low to be recognized by the test. In case of suspicion, repeat the test after 1-2 days, as the SARS-CoV-2/Influenza virus/RSV/Adenovirus cannot be precisely detected in all phases of an infection.

6. If the test result is positive, what should I do?

A positive result indicates that you (or the patient) probably have (has) a SARS-CoV-2, Influenza A, Influenza B, RSV or Adenovirus infection. Please see a doctor for medical aid.

PERFORMANCE CHARACTERISTICS

Accuracy

The test was evaluated with clinical specimens. RT-PCR was used as the reference method. Study results are presented in the tables below.

COVID-19 Test:				
Method	RT-PCR			
	Positive	Negative	Total	
Rapid Test	Positive	213	1	214
	Negative	5	1377	1382
Total	218		1378	1596
Sensitivity	97.7% (95% CI*: 94.7%~99.3%)			
Specificity	99.9% (95% CI*: 99.6%~100.0%)			
Accuracy	99.6% (95% CI*: 99.2%~99.9%)			

FLU A+B Test:

Method	FLU A			FLU B		
	RT-PCR			RT-PCR		
Rapid Test	Positive	Negative	Total	Positive	Negative	Total
	179	1	180	137	0	137
Negative		939	948	8	983	991

Total	188	940	1128	145	983	1128
Sensitivity	95.2% (95% CI*: 91.1%~97.8%)			94.5% (95% CI*: 89.4%~97.6%)		
Specificity	99.9% (95% CI*: 99.4%~100.0%)			>99.9% (95% CI*: 99.6%~100.0%)		
Accuracy	99.1% (95% CI*: 98.4%~99.6%)			99.3% (95% CI*: 98.6%~99.7%)		

RSV Test:

Method		RT-PCR		Total
		Positive	Negative	
Rapid Test	Positive	151	2	153
	Negative	6	833	839
Total		157	835	992
Sensitivity		96.2% (95% CI*: 91.9%~98.6%)		
Specificity		99.8% (95% CI*: 99.1%~100.0%)		
Accuracy		99.2% (95% CI*: 98.4%~99.7%)		

Adenovirus Test:

Method		RT-PCR		Total
		Positive	Negative	
Rapid Test	Positive	118	2	120
	Negative	3	789	792
Total		121	791	912
Sensitivity		97.5% (95% CI*: 92.9%–99.5%)		
Specificity		99.8% (95% CI*: 99.1%–100.0%)		
Accuracy		99.5% (95% CI*: 98.7%–99.8%)		

* Confidence Intervals

Sensitivity

The COVID-19/FLU A+B/RSV/Adenovirus Antigen Rapid Test can detect various virus strains of SARS-CoV-2, Influenza A/B viruses, RSV and Adenovirus. The confirmed LODs for each virus strain are presented in the table below.

Virus Strain	LOD
WHO SARS-CoV-2 antigen International Standard	30 IU/mL
BetaCoV/Wuhan/IPBCAMS-WH-01/2019	78 TCID ₅₀ /mL
A/Sydney/5/2021(H1N1)	50 TCID ₅₀ /mL
A/California/07/2009(H1N1)	5.25 × 10 ³ CEID ₅₀ /mL
A/South Australia/69/2019(H3N2)	50 TCID ₅₀ /mL
A/Wisconsin/15/2009(H3N2)	6.2 × 10 ⁴ CEID ₅₀ /mL
A/Anhui/1/2013(H7N9)	650 TCID ₅₀ /mL
A/Guangdong/175F003/2016(H7N9)	4.5 × 10 ⁴ TCID ₅₀ /mL
B/Austria/1359417/2021(Victoria)	50 TCID ₅₀ /mL
B/Florida/78/2015(Victoria)	1.4 × 10 ⁴ TCID ₅₀ /mL
B/Darwin/58/2019(Yamagata)	100 TCID ₅₀ /mL
B/Brisbane/1/2007(Yamagata)	3.16 × 10 ⁴ TCID ₅₀ /mL
A2(RSV type A)	2.5 × 10 ⁴ TCID ₅₀ /mL
B/VW/14617/85(RSV type B)	1.0 × 10 ⁴ TCID ₅₀ /mL
GB(Adenovirus type 3)	1.0 × 10 ⁴ TCID ₅₀ /mL
Gömen(Adenovirus type 7)	5.0 × 10 ⁴ TCID ₅₀ /mL

The test can also detect strains of SARS-CoV-2, Influenza A/B viruses, RSV and Adenovirus as shown in the following table.

Virus Strains	Subtype	Concentration (TCID ₅₀ /mL)
SARS-CoV-2	Delta variants (B.1.617.2)	1.0 × 10 ⁴
SARS-CoV-2	Delta variants (B.1.617.3)	1.0 × 10 ⁴
SARS-CoV-2	Delta variants (AY Sub-lineages)	1.0 × 10 ⁴
SARS-CoV-2	Omicron variants (B.1.1.529)	1.0 × 10 ⁴
SARS-CoV-2	Omicron variants (BA.2)	1.0 × 10 ⁴
A/Brisbane/02/2018	Influenza A(H1N1)	1.5 × 10 ⁴
A/Victoria/2570/2019	Influenza A(H1N1)	2.2 × 10 ⁴
A/Sydney/175/2022	Influenza A(H1N1)	2.9 × 10 ⁴
A/Darwin/122/2018	Influenza A(H1N1)	3.3 × 10 ⁴
A/Darwin/6/2018	Influenza A(H1N1)	1.1 × 10 ⁴
A/Hangzhou/A158/2023	Influenza A(O9H1N1)	5.25 × 10 ⁶
A/Hangzhou/A419/2024	Influenza A(O9H1N1)	2.63 × 10 ⁶
A/Hubei/WL1481/2023	Influenza A(O9H1N1)	1.51 × 10 ⁶
A/Brisbane/192/2017	Influenza A(H3N2)	3.1 × 10 ⁴
A/Perth/9/2019	Influenza A(H3N2)	1.9 × 10 ⁴
A/Darwin/9/2021	Influenza A(H3N2)	4.4 × 10 ⁴
A/Victoria/6/2022	Influenza A(H3N2)	3.7 × 10 ⁴
A/Sydney/1020/2018	Influenza A(H3N2)	6.4 × 10 ⁴
A/Hubei/563/2019	Influenza A(H3N2)	5.9 × 10 ⁴
A/Hangzhou/B870/2023	Influenza A(H3N2)	1.47 × 10 ⁷
A/Hubei/VQ180801/2019	Influenza A(H5N1)	2.5 × 10 ⁷
A/Anhui/VQ190501/2019	Influenza A(H5N1)	2.0 × 10 ⁷
A/Hunan/VQ180901/2019	Influenza A(H5N1)	3.1 × 10 ⁷
A/Hubei/c18259255/2018	Influenza A(H5N6)	1.0 × 10 ⁵
A/LVS/2/2018	Influenza A(H5N6)	1.2 × 10 ⁶
A/Hubei/c18255022/2018	Influenza A(H7N9)	1.0 × 10 ⁴
A/Shanghai/02/2013	Influenza A(H7N9)	1.5 × 10 ⁴
A/Hubei/c180510359/2018 (H9N2)	Influenza A(H9N2)	1.0 × 10 ⁴
A/LVS/2/2018	Influenza A (H9N2)	1.4 × 10 ⁶
B/Victoria/2110/19	Influenza B(Victoria)	1.7 × 10 ⁴
B/Brisbane/35/2018	Influenza B(Victoria)	1.5 × 10 ⁴
B/Victoria/705/2018	Influenza B(Victoria)	5.7 × 10 ⁴
B/Hangzhou/B269/2024	Influenza B(Victoria)	6.67 × 10 ⁶
B/Hubei/1832/2016	Influenza B(Victoria)	4.74 × 10 ⁶
B/South Australia/67/2018	Influenza B(Yamagata)	2.4 × 10 ⁴
B/Victoria/706/2018	Influenza B(Yamagata)	8.3 × 10 ⁶
B/Hangzhou/B52/2018	Influenza B(Yamagata)	5.62 × 10 ⁶

B/Hubei/24/2018	Influenza B(Yamagata)	2.16 × 10 ⁶
Long	RSV type A	1.6 × 10 ⁴
HRSV-A/Hep2/Hubei/2024-XN-335/2024	RSV type A	1.1 × 10 ⁴
18537	RSV type B	1.6 × 10 ⁴
RSVB/Human/Hubei/201710052/2018	RSV type B	1.0 × 10 ⁴
Adenoid 71	Human adenovirus 1	1.0 × 10 ⁷
Adenoid 6	Human adenovirus 2	1.0 × 10 ⁶
RI-67	Human adenovirus 4	5 × 10 ⁶
Clinical Virus Isolate	Human adenovirus 5	4 × 10 ⁶
Tonsil 99	Human adenovirus 6	5 × 10 ⁷
7A	Human adenovirus 7	3.16 × 10 ³
Clinical Virus Isolate	Human adenovirus 11	3.1 × 10 ⁴
Clinical Virus Isolate	Human adenovirus 14	3.1 × 10 ⁴
Clinical Virus Isolate	Human adenovirus 55	1.6 × 10 ⁴

Hook Effect

SARS-CoV-2 positive specimen of 1.0 × 10⁷ TCID₅₀/mL, Influenza A virus positive specimen of 4.22 × 10⁷ TCID₅₀/mL, Influenza B virus positive specimen of 1.78 × 10⁷ TCID₅₀/mL, RSV positive specimen of 1.6 × 10⁷ TCID₅₀/mL and Adenovirus positive specimen of 1.6 × 10⁴ TCID₅₀/mL were tested and no high-dose hook effect was observed.

Precision-Repeatability

Precision-repeatability was determined by using specimens: negative, SARS-CoV-2 weak positive, SARS-CoV-2 strong positive, Influenza A weak positive, Influenza A strong positive, Influenza B weak positive, Influenza B strong positive, RSV weak positive, RSV strong positive, Adenovirus weak positive and Adenovirus strong positive specimens. The study was performed 5 replicates per day for 5 consecutive days by one operator using 1 lot of the test. No difference was detected in intra lot.

Precision-Reproducibility

Precision-reproducibility was determined using specimens: negative, SARS-CoV-2 weak positive, SARS-CoV-2 strong positive, Influenza A weak positive, Influenza A strong positive, Influenza B weak positive, Influenza B strong positive, RSV weak positive, RSV strong positive, Adenovirus weak positive and Adenovirus strong positive specimens. The study was performed 5 replicates per day for 5 consecutive days in 3 different sites using

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EN Streptococcus A

INTENDED USE

The test is a rapid chromatographic immunoassay for the qualitative detection of Streptococcal A antigens from self-collected throat swab specimens to aid in the diagnosis of Group A Streptococcal infection. The test is not automated. For self-testing use. For *in vitro* diagnostic use only.

HOW TO USE THIS TEST

- The **ZuhauseTEST Streptococcus A** is an immunological test designed to specifically identify Streptococcus A antigens in pharyngeal infections. It helps to understand quickly if a sore throat is caused by Group A Streptococcus.
- If the test result is negative but you continue to have symptoms, you may not have Strep A infection. However you should ask your healthcare provider for advice.
- If your test result is positive, it means antigens from the bacterium that causes Strep A infection have been found in your specimen, and you are likely to have a Strep A infection.

Zuhause

Rapid test for the qualitative detection of Streptococcal A antigens in throat swabs

FR Streptocoque A

Instructions for use

Group B Streptococcus	Group C Streptococcus
Group F Streptococcus	Group G Streptococcus
<i>Streptococcus pneumoniae</i>	<i>Streptococcus sanguinis</i>
<i>Streptococcus mutans</i>	<i>Staphylococcus epidermidis</i>
<i>Staphylococcus aureus</i>	<i>Serratia marcescens</i>
<i>Corynebacterium diphtheriae</i>	<i>Klebsiella pneumoniae</i>
<i>Candida albicans</i>	<i>Bordetella pertussis</i>
<i>Enterococcus faecalis</i>	<i>Neisseria gonorrhoeae</i>
<i>Neisseria meningitidis</i>	<i>Neisseria subflava</i>
<i>Neisseria sicca</i>	<i>Haemophilus influenzae</i>
<i>Branhamella catarrhalis</i>	<i>Pseudomonas aeruginosa</i>

Interfering Substances: The following compounds were tested using the **ZuhauseTEST Streptococcus A** and no interference was observed.

Substance	Concentration
Cherry Halls cough drops	5 mg/mL
Menthhol Halls cough drops	5 mg/mL
Robitussin cough syrup	10% (vol/vol)
Dimetapp cough syrup	10% (vol/vol)
Vicks Chloraseptic spray	20% (vol/vol)
Cepacol Chloraseptic spray	20% (vol/vol)
Listerine mouthwash	20% (vol/vol)
Scope mouthwash	20% (vol/vol)
Blood (human)	4% (vol/vol)
Mucin	0.5%
Saliva	5%

STORAGE AND STABILITY

- Store the kit at 2-30°C. DO NOT FREEZE.
- Keep the kit away from sunlight, moisture and heat.
- Do not use the kit beyond the expiration date.
- Please use the test cassette within one hour after removing it from the foil pouch at room temperature (15-30°C).

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- For self-testing *in vitro* diagnostic use only. Do not use it after the expiration date.
- Store the kit in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. Do not use it if any of the test kit content or packaging is damaged or opened.
- Use the test only once and follow the test procedures strictly. Do not dismantle the test cassette or touch the test window of the test cassette.
- Keep the kit out of the reach of children. Test for children should be done by an adult.
- For people with reading difficulties, please use the kit with the assistance of another adult.
- The Buffer 1 contains 13.8% NaNO₂. It is harmful if swallowed. IF SWALLOWED, rinse the mouth. Call a POISON CENTER doctor if you feel unwell.
- Do not swallow Buffer 1, Buffer 2 or any other components in the box!
- If the buffer solutions contact the skin or eye, flush it with large volume of water.
- Buffer 1 should be stored away from light and be kept in the black packaging bag before use.
- Wash hands thoroughly before and after handling.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- The used kit should be discarded according to local regulations.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.
- Components provided in the kit are approved for use in the **ZuhauseTEST Streptococcus A**. Do not use any component in other commercial kits.
- This test provides only a preliminary test result and does not replace medical advice. Always seek guidance from a qualified healthcare provider before making any health-related decisions.

LIMITATIONS

- This test only indicates the presence of Group A Streptococcus bacteria in throat swab specimens. Neither the quantitative value nor the rate of increase in Streptococcus A antigen concentration can be determined by this qualitative test.
- A false negative result may be obtained if the concentration of the Streptococcus A antigen present in the throat swab is not adequate or below the detectable level of the test.
- Too much blood or mucus on the swab specimen may interfere with test performance and may yield a false positive result. Avoid touching the tongue, cheeks, teeth, or any bleeding areas of the mouth with the swab when collecting specimens.
- The **ZuhauseTEST Streptococcus A** provides preliminary test results and should not be used as the sole basis for treatment or other management decision. As with all diagnostic tests, a confirmed diagnosis should only be made after evaluating other clinical information available to the physician.
- A positive result for Streptococcus A does not preclude an underlying co-infection with another pathogen; therefore, the possibility of an underlying bacterial infection should be considered.
- False negative or false positive results may occur if a specimen is improperly collected or handled.

BIBLIOGRAPHY

- Hammad S, Kawabata S, Nakagawa I (2015) Molecular and genomic characterization of pathogenic traits of group A *Streptococcus pyogenes*. Proc Jpn Acad Ser B, 91(10) 539–559.
- Webb K H (1998) Does culture confirmation of high-sensitivity rapid streptococcal tests make sense? A medical decision analysis. Pediatrics, 101(2) e2.
- Bisno A L, Gerber M A, Kaplan E L, Schwartz R H (1997) Diagnosis and management of group A streptococcal pharyngitis: A practice guideline. Clin Infect Dis 25(3) 574–583.

QUESTION & ANSWER

- How does the ZuhauseTEST Streptococcus A work?**
The test can qualitatively detect Streptococcus A antigens in self-collected swab specimens through the specific antibodies it contains. A positive result indicates that Streptococcus A antigens are probably present in the specimen.

- When should the test kit be used?**
You can perform this test at home as an early diagnosis of group A streptococcal infection if you have the symptoms such as sore throat, red and swollen tonsils, small red patches on the back of the palate (soft or hard), swollen lymph nodes, fever, headache, nausea or vomiting, especially in children.
- Can the result be incorrect?**
The result is accurate, as long as the instructions are carefully followed. Nevertheless, the result can be incorrect if inadequate or improper specimen is collected or the **ZuhauseTEST Streptococcus A** gets wet before testing, or if the number of drops added is less than 3 or more than 4. Besides, due to immunological principles involved, there exist the chances of false positive or false negative results in rare cases under proper use and normal conditions. A consultation with the doctor is always recommended for such tests based on immunological principles.
- How do I interpret the test result if the color and the intensity of the lines are different?**
The color and the intensity of the lines have no importance for test result interpretation. The lines should only be homogeneous and clearly visible.
- My test result is negative. Does that mean I'm not infected?**
A negative result means that you are probably not infected or that the amount of bacteria is too small to be recognized by the test. In case of suspicion, repeat the test after 1-2 days, as the Strep A cannot be precisely detected in all phases of an infection.
- If the test result is positive, what should I do?**
A positive result indicates that you may have group A streptococcal infection and please see a doctor for medical aid. Please be aware that the **ZuhauseTEST Streptococcus A** cannot distinguish between live bacteria and dead bacteria. People who have recently recovered from Streptococcus A infection may test positive because of inactive bacteria that are still present in the body.

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EN A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

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Steriler Tupfer/Sterile swab/écouvillon stérile/hisopo estéril:

Hangzhou Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangting District,
Yangzhou, Jiangsu 225109 China

oder/ou/ou/o

Jiangsu Hanheng Medical Technology Co., Ltd
#118 Xinyuan 3rd Road, Nxbel District, Changzhou,
213031, Jiangsu P.R. China

Zungenspatel/Tongue depressor/
abaisse-langue/depresor de lengua:

Dalian Rongbang Medical Healthy Devices Co.,Ltd
Maoyingzi Hamlet, Dalianwan Town,
Ganjingzi District, Dalian, 116113

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EN If only one colored line appears in the control line region (C) and no colored line appears in the test line region (T), the result is negative. It indicates that no Streptococcus A antigen has been detected in the specimen or that the amount is too small to be detected.

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- Hammada S, Kawabata S, Nakagawa I (2015) Molecular and genomic characterization of pathogenic traits of group A Streptococcus pyogenes. Proc Jpn Acad Ser B, 91(10) 539-559.
- Webb K H (1998) Does culture confirmation of high-sensitivity rapid streptococcal tests make sense? A medical decision analysis. Pediatrics, 101(2) e2.
- Bisno A L, Gerber M A, Kaplan E L, Schwartz R H (1997) Diagnosis and management of group A streptococcal pharyngitis: A practice guideline. Clin Infect Dis 25(3) 574-583.

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EN If the line in the control line region (C) fails to appear, the result is invalid. Incorrect testing procedures or test reagent deterioration are the most likely reasons. Please review the testing procedures and repeat the test with a new test kit. If the problem persists, discontinue using the kit immediately and contact your local distributor.

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

EN Allow the test, specimen and buffer to reach room temperature (15-30°C) prior to testing. Please wash and dry your hands thoroughly before testing. Read the package insert carefully.

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
EN Withdraw the swab and be careful not to touch the tongue, cheeks and teeth.

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
EN Start a timer and read the result 5 minutes later after adding the specimen. Do not interpret the result after 10 minutes.

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
EN Remove the cover of the Buffer 2 tube and insert the tube into the dotted circle on the packaging box. Break off the Buffer 1 at the breaking point (the narrowest place). Then add entire Buffer 1 solution to the Buffer 2 tube. The mixed solution should be used immediately.

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
EN Place the swab immediately into the tube containing mixed solution of Buffer 1 and Buffer 2. Then rotate the swab against the inside of the tube about 15 times. Leave the swab in the tube for 1 minute.

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
EN When finished, dispose of all used components according to local regulations.

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EN Open the package containing the sterile swab. Remove the swab and avoid touching the soft tip.

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Pantone 7461C

CMYK

Attention: By approving the enclosed design draft, you (**AllTest's** Customer) accept all responsibility for the accuracy of the design. If an error is detected following the printing or manufacturing of a material, you (**AllTest's** Customer) are responsible for the cost of any inventory which is deemed unsuitable for sale and /or not compliant with local regulations.

Part Number	14603431300	Description	(彩色)Zuhause TEST ISTB-N502H(630000) CE0123 4国语言说明书(德英法西)(欧代新名称)(For A1703)	
Size	420x297mm	Material	70g双胶纸,双面,折好到货,折后尺寸约53x73mm	
Designer	zhu	Design Date/Version	Apr 30 2026/K	Mold Num.
Confirmation and Signature of Customer/Date				