Zuhause TEST

DE Streptokokken A

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Zuhause TEST

EN Streptococcus A

Rapid test for the qualitative detection of Streptococcal A antigens in throat swabs

INSTRUCTIONS FOR USE

HOW IT WORKS

Streptococcus pyogenes (S. pyogenes) is gram-positive, nonmotile, facultatively anaerobic cocci. It is a human adapted β-hemolytic, Lancefield group A Streptococcus, commonly referred to as GAS. GAS usually colonizes the throat or skin epithelial surfaces and causes a wide variety of clinical manifestations such as noninvasive pharyngitis, dermatitis, and scarlet fever, as well as invasive systemic infections such as necrotizing fasciitis (NF) and streptococcal toxic shock syndrome (STSS) in humans.[1] Left untreated, these infections can lead to serious complications, including rheumatic fever and peritonsillar abscess.[2] Traditional identification procedures for Group A Streptococci infection involve the isolation and identification of viable organisms using techniques that require 24 to 48 hours or longer.[3]
The **ZuhauseTEST Streptococcus A** is a qualitative, lateral flow immunoassay for the detection of Strep A antigens in throat swab specimens. When the extracted specimen is added into the specimen well, Strep A antigens, if present in the specimen, reacts with an antibody to Strep A that is coated onto latex particles, forming a mixture. The mixture migrates up the membrane and reacts with the antibody to Strep A on the membrane, and generate a red line in the test line region (T) as a result. The presence of this red line in the test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a blue line will always appear in the control line region (C), indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test contains Anti-Streptococcus Group A antibody, and Streptavidin-Rabbit IgG coated on the membrane; and contains Anti-Streptococcus Group A antibody coated particles and Biotin-BSA coated particles.

MATERIALS

- 1 Test cassette [Strep A]
- 1 Buffer 1 [0.24 mL (13.8% $NaNO_2$, 0.0004% Phenol red)]
- 1 Buffer 2 [0.16 mL (0.5184% Citric acid, 0.02% Proclin 300)]
- 1 Sterile swab
- 1 Tongue depressor
- 1 Instruction for use

Materials Required but Not Provided:

- 1 timer

PERFORMANCE CHARACTERISTICS

Accuracy: A study has been conducted on altogether 304 throat swab specimens. The **ZuhauseTEST Streptococcus A** was studied in parallel comparison with Culture.

Method Comparison		Professional		
		Culture		Total
Zuhause-TEST Streptococcus A	Positive	Positive	Negative	
		96	4	100
	Negative	5	199	204
Total		101	203	304

Sensitivity: 95.1% (95% CI*: 88.8%-98.4%)
Specificity: 98.0% (95% CI*: 95.0%-99.5%)
Accuracy: 97.0% (95% CI*: 94.5%-98.6%)
*Confidence Interval

Lay User Clinical Study: A lay user clinical study was conducted to evaluate the clinical performance and usability of **ZuhauseTEST Streptococcus A** in the hands of lay users. A total of 100 lay users were included in this study, including 31 positive cases and 69 negative cases. The study results indicated that the positive percent agreement, the negative percent agreement and the overall percent agreement of **ZuhauseTEST Streptococcus A** are 96.8%, 100.0%, and 99.0% respectively when comparing with another rapid test, demonstrating favorable consistency between this product and the reference kit. In addition, the study results also proved that the instructions in package Insert, labeling are clear and easy to understand, and the rapid test is easy to use by lay users. Thus, this **ZuhauseTEST Streptococcus A** is safe, suitable and effective for lay users in detecting Strep A antigens in throat swab specimens in a home setting.

Method Comparison		Another Rapid Test		Total
Zuhause-TEST Streptococcus A	Positive	Positive	Negative	
		30	0	30
	Negative	1	69	70
Total		31	69	100

Positive Percent Agreement:
96.8% (95%CI*: 83.8%-99.4%)
Negative Percent Agreement:
100.0% (95%CI*: 94.7%-100.0%)
Overall Percent Agreement:
99.0% (95%CI*: 94.6%-99.8%)
*Confidence Interval

Analytical Sensitivity: The **ZuhauseTEST Streptococcus A** can detect Streptococcus A of as low as $1.0X10^7$ CFU/mL (10^5 CFU/test).
Hook Effect: There is no dose hook effect with the **ZuhauseTEST Streptococcus A** when testing Streptococcus A of no more than $1.0X10^{12}$ CFU/mL ($1.0X10^{10}$ CFU/test).
Precision-Repeatability: Repeatability was determined using four specimens: a negative specimen (0.5% BSA-PBS), and three Streptococcus A positive specimens ($1.0X10^7$ CFU/mL, $2.5X10^7$ CFU/mL, $1.0X10^8$ CFU/mL). The study was performed by testing each specimen in 5 replicates per day for 5 consecutive days by one operator using 1 lot of **ZuhauseTEST Streptococcus A**. No difference was detected in intra lot.
Precision-Reproducibility: Reproducibility was determined using seven specimens: a negative specimen (0.5% BSA-PBS), and six Streptococcus A positive specimens ($5.0X10^6$ CFU/mL, $1.0X10^7$ CFU/mL, $1.5X10^7$ CFU/mL, $2.0X10^7$ CFU/mL, $2.5X10^7$ CFU/mL, $1.0X10^8$ CFU/mL). The study was performed by testing each specimen in 6 replicates per day for 5 consecutive days in 3 different sites using 3 separate lots of **ZuhauseTEST Streptococcus A** by three operators (one lot/operator per site). The study results indicated high precision of this rapid test.
Cross Reactivity: The following organisms were tested at $1.0X10^7$ CFU/mL with the **ZuhauseTEST Streptococcus A** and were all found to be negative. No mucoid-producing strains were tested.

Group B Streptococcus	Group C Streptococcus
Group F Streptococcus	Group G Streptococcus
Streptococcus pneumoniae	*Streptococcus sanguinis*
Streptococcus mutans	*Staphylococcus epidermidis*
Staphylococcus aureus	*Serratia marcescens*
Corynebacterium diphtheriae	*Klebsiella pneumoniae*
Candida albicans	*Bordetella pertussis*
Enterococcus faecalis	*Neisseria gonorrhoeae*
Neisseria meningitidis	*Neisseria subflava*
Neisseria sicca	*Haemophilus influenzae*
Branhamella catarrhalis	*Pseudomonas aeruginosa*

Interfering Substances: The following compounds were tested using the **ZuhauseTEST Streptococcus A** and no interference was observed.

Substance	Concentration
Cherry Halls cough drops	5 mg/mL
Menthol Halls cough drops	5 mg/mL
Robitussin cough syrup	10% (vol/vol)
Dimetapp cough syrup	10% (vol/vol)
Vicks Chloraseptic spray	20% (vol/vol)
Cepacol Chloraseptic spray	20% (vol/vol)
Listerine mouthwash	20% (vol/vol)
Scope mouthwash	20% (vol/vol)
Blood (human)	4% (vol/vol)
Mucin	0.5%
Saliva	5%

STORAGE AND STABILITY

1. Store the kit at 2-30°C. DO NOT FREEZE.
2. Keep the kit away from sunlight, moisture and heat.
3. Do not use the kit beyond the expiration date.
4. Please use the test cassette within one hour after removing it from the foil pouch at room temperature (15-30°C).

WARNINGS AND PRECAUTIONS

Please read all the information in this package insert before performing the test.

1. For self-testing *in vitro* diagnostic use only. Do not use it after the expiration date.
2. Store the kit in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. Do not use it if any of the test kit content or packaging is damaged or opened.
3. Use the test only once and follow the test procedures strictly. Do not dismantle the test cassette or touch the test window of the test cassette.
4. Keep the kit out of the reach of children. Test for children should be done by an adult.
5. For people with reading difficulties, please use the kit with the assistance of another adult.
6. The Buffer 1 contains 13.8% $NaNO_2$. It is harmful if swallowed. IF SWALLOWED, rinse the mouth. Call a POISON CENTER doctor if you feel unwell.
7. Do not swallow Buffer 1, Buffer 2 or any other components in the box!
8. If the buffer solutions contact the skin or eye, flush it with large volume of water.
9. Buffer 1 should be stored away from light and be kept in the black packaging bag before use.
10. Wash hands thoroughly before and after handling.
11. Do not eat, drink or smoke in the area where the specimens or kits are handled.
12. The used kit should be discarded according to local regulations.
13. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.
14. Components provided in the kit are approved for use in the **ZuhauseTEST Streptococcus A**. Do not use any component in other commercial kits.
15. This test provides only a preliminary test result and does not replace medical advice. Always seek guidance from a qualified healthcare provider before making any health-related decisions.

LIMITATIONS

1. This test only indicates the presence of Group A Streptococcus bacteria in throat swab specimens. Neither the quantitative value nor the rate of increase in Streptococcus A antigen concentration can be determined by this qualitative test.
2. A false negative result may be obtained if the concentration of the Streptococcus A antigen present in the throat swab is not adequate or below the detectable level of the test.
3. Too much blood or mucus on the swab specimen may interfere with test performance and may yield a false positive result. Avoid touching the tongue, cheeks, teeth, or any bleeding areas of the mouth with the swab when collecting specimens.
4. The **ZuhauseTEST Streptococcus A** provides preliminary test results and should not be used as the sole basis for treatment or other management decision. As with all diagnostic tests, a confirmed diagnosis should only be made after evaluating other clinical information available to the physician.
5. A positive result for Streptococcus A does not preclude an underlying co-infection with another pathogen; therefore, the possibility of an underlying bacterial infection should be considered.
6. False negative or false positive results may occur if a specimen is improperly collected or handled.

BIBLIOGRAPHY

1. Hammada S, Kawabata S, Nakagawa I (2015) Molecular and genomic characterization of pathogenic traits of group A *Streptococcus pyogenes*. Proc Jpn Acad Ser B, 91(10) 539–559.
2. Webb K H (1998) Does culture confirmation of high-sensitivity rapid streptococcal tests make sense? A medical decision analysis. Pediatrics, 101(2) e2.
3. Bisno A L, Gerber M A, Kaplan E L, Schwartz R H (1997) Diagnosis and management of group A streptococcal pharyngitis: A practice guideline. Clin Infect Dis 25(3) 574–583.

QUESTION & ANSWER

1. **How does the ZuhauseTEST Streptococcus A work?**
 The test can qualitatively detect Streptococcus A antigens in self-collected swab specimens through the specific antibodies it contains. A positive result indicates that Streptococcus A antigens are probably present in the specimen.
2. **When should the test kit be used?**
 You can perform this test at home as an early diagnosis of group A streptococcal infection if you have the symptoms such as sore throat, red and swollen tonsils, small red patches on the back of the palate (soft or hard), swollen lymph nodes, fever, headache, nausea or vomiting, especially in children.
3. **Can the result be incorrect?**
 The result is accurate, as long as the instructions are carefully followed. Nevertheless, the result can be incorrect if inadequate or improper specimen is collected or the **ZuhauseTEST Streptococcus A** gets wet before testing, or if the number of drops added is less than 3 or more than 4. Besides, due to immunological principles involved, there exist the chances of false positive or false negative results in rare cases under proper use and normal conditions. A consultation with the doctor is always recommended for such tests based on immunological principles.
4. **How do I interpret the test result if the color and the intensity of the lines are different?**
 The color and the intensity of the lines have no importance for result interpretation. The lines should only be homogeneous and clearly visible.
5. **My test result is negative. Does that mean I'm not infected?**
 A negative result means that you are probably not infected or that the amount of bacteria is too small to be recognized by the test.
 In case of suspicion, repeat the test after 1-2 days, as the Strep A cannot be precisely detected in all phases of an infection.
6. **If the test result is positive, what should I do?**
 A positive result indicates that you may have group A streptococcal infection and please see a doctor for medical aid.
 Please be aware that the **ZuhauseTEST Streptococcus A** cannot distinguish between live bacteria and dead bacteria. People who have recently recovered from Streptococcus A infection may test positive because of inactive bacteria that are still present in the body.

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EN A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

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Steriler Tupfer/Sterile swab/écouvillon stérile/hisopo estéril:
Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,Jiangsu 225109 China
CE 0197
oder/or/ou/o
Jiangsu Hanheng Medical Technology Co., Ltd
#118 Xinyuan 3rd Road, Xnbei District, Changzhou, 213031, Jiangsu P.R. China
CE 0197
Zungenspatel/Tongue depressor/abaisse-langue/depresor de lengua:
Dalian Rongbang Medical Healthy Devices Co.,Ltd
Maoyingzi Hamlet, Dalianwan Town, Ganjingzi District, Dalian, 116113
CE 0197

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EN If two colored lines appear both in the control line region (C) and the test line region (T), the result is positive. It indicates that you may be in a stage of Streptococcus A infection and should consult with your doctor.
***NOTE:** The line in the control line region (C) is blue and the line in the test line region (T) is red. In addition, the color intensity of the line in the test line region will vary based on the amount of Strep A antigen present in the specimen. Thus any shade of a colored line in test line region indicates a positive result.

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EN If only one colored line appears in the control line region (C) and no colored line appears in the test line region (T), the result is negative. It indicates that no Streptococcus A antigen has been detected in the specimen or that the amount is too small to be detected.

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EN If the line in the control line region (C) fails to appear, the result is invalid. Incorrect testing procedures or test reagent deterioration are the most likely reasons. Please review the testing procedures and repeat the test with a new test kit. If the problem persists, discontinue using the kit immediately and contact your local distributor.

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

EN Allow the test, specimen and buffer to reach room temperature (15-30°C) prior to testing. Please wash and dry your hands thoroughly before testing. Read the package insert carefully.
Note: Check kit components and confirm the kit is not expired prior to use.

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
EN Remove the cover of the Buffer 2 tube and insert the tube into the dotted circle on the packaging box. Break off the Buffer 1 at the breaking point (the narrowest place). Then add entire Buffer 1 solution to the Buffer 2 tube. The mixed solution should be used immediately.
Note: Buffer 1 needs to stay in the black packaging bag before use.

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
EN Open the package containing the sterile swab. Remove the swab and avoid touching the soft tip.

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
EN Stand in front of a mirror, tilt your head backwards slightly and open your mouth as wide as possible. Depress the tongue and hold it down with a tongue depressor or a spoon. Collect the specimen by rubbing the sterile swab tip on the surface of one or both tonsils or the posterior pharyngeal wall.
Note: Other areas of the oral pharynx and mouth are not acceptable specimen collection sites, and could lead to false negative results.

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
EN Withdraw the swab and be careful not to touch the tongue, cheeks and teeth.
Note: Testing should be performed immediately after the specimen has been collected.

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
EN Place the swab immediately into the tube containing mixed solution of Buffer 1 and Buffer 2. Then rotate the swab against the inside of the tube about 15 times. Leave the swab in the tube for 1 minute.

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
EN Remove the swab while squeezing the swab head against the inside of the tube to release as much liquid as possible from the swab.
Note: Discard the swab in compliance with local regulations.

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EN Fit the tube tip onto the tube. Take the test cassette out from the foil pouch and place it on a flat and level surface. Invert the tube and add 3 drops of solution into the specimen well (S) of the test cassette.
Note: If a drop of solution added contains air bubbles, add another drop to the well.

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
EN Start a timer and read the result 5 minutes later after adding the specimen. Do not interpret the result after 10 minutes.
Note: Do not touch or move the test cassette during the waiting time.

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EN When finished, dispose of all used components according to local regulations.

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